

## T-Pro Transfection Reagent II



Store at  
2~8°C

(JT98-T002S) 0.5 ml  
(JT98-T002M) 1.0 ml

**This product is for laboratory research ONLY and not for diagnostic use.**

<b>Description</b>	T-Pro Transfection Reagent II is a proprietary formulation for the transfection of DNA and RNA into eukaryotic cells and forms the covalent binding with nucleic acid to provide the highly efficiency and low cyto-toxicity. T-Pro Transfection Reagent II can be used in many cell types including HEK293, 293T, 293E, CHO, COS1, HeLa, NIH 3T3, Sf9 and Sf21. Furthermore, T-Pro Transfection Reagent II also has high transfection efficiency when applied in other eukaryotic cells
<b>Properties</b>	<ul style="list-style-type: none"><li>- Low cytotoxicity for most of tumor cell lines and primary cells.</li><li>- Exceptional transfection efficiency of a broad range of cell types.</li><li>- Efficient transfection with or without serum.</li><li>- High levels of recombinant protein production.</li><li>- Simple, robust transfection procedure.</li></ul>
<b>Storage</b>	T-Pro Transfection Reagent II is stable for 2~8°C

### Procedure

#### For adhesion cells

Note: Use the following procedure to transfect mammalian cells in a 6-well plate format. For other formats, see Scaling Up or Down Transfections. All amounts and volumes are given on a per well basis.

- 1 One day before transfection, plate  $0.5-2 \times 10^5$  cells in 500  $\mu$ l of growth medium without antibiotics so that cells will be 70-95% confluent at the time of transfection (table 1).
- 2 For each transfection sample, prepare complexes as follows T-Pro Transfection Reagent II ( $\mu$ l) : DNA ( $\mu$ g) (3 : 1~12 : 1)
  - a. Dilute DNA 2.5 $\mu$ l (1 $\mu$ g/ $\mu$ l) in 25  $\mu$ l of serum free medium or HBS buffer (pH 7.4, 150 mM NaCl, 20mM HEPES) Mix gently.
  - b. Mix T-Pro Transfection Reagent III 15 $\mu$ l gently before use, and then dilute the appropriate amount in 50  $\mu$ l of serum free medium or HBS buffer (pH 7.4, 150 mM NaCl, 20mM HEPES) Mix gently. Incubate for 5 minutes at room temperature.  
**Note: Combine diluted T-Pro Transfection Reagent II with diluted DNA within 30 minutes.**
  - c. After 5 minute incubation, combine the diluted DNA with diluted T-Pro Transfection Reagent II (total volume = 100  $\mu$ l). Mix gently and incubate for 15 minutes at room temperature (solution may appear cloudy).
- 3 Add the 100  $\mu$ l of complexes to each well containing cells and medium. Mix gently by rocking the plate back and forth.
- 4 Incubate cells at 37°C in a CO<sub>2</sub> incubator for 24-48 hours prior to testing for transgene expression. It is not necessary to change the medium, but medium may be replaced after 6-24 hours.

**Table 1.** A Guideline for Seeding Adherent Cells Prior to Transfection in Different Culture Formats.

Culture Dishes	Surface Area (cm <sup>2</sup> )	Number of Cells to Seed
T175 Flask	175	0.7 – 1.4 x 10 <sup>7</sup>
T75 Flask	75	3.0 – 6.0 x 10 <sup>6</sup>
100 mm Dish	58	2.2 – 4.4 x 10 <sup>6</sup>
60 mm Dish	21	0.9 – 1.8 x 10 <sup>6</sup>
35 mm Dish	9.6	3.5 – 7.0 x 10 <sup>5</sup>
6-well Plate	9.6	4.0 – 8.0 x 10 <sup>5</sup>
12-well Plate	3.5	1.5 – 3.0 x 10 <sup>5</sup>
24-well Plate	1.9	0.8 – 1.6 x 10 <sup>5</sup>
48-well Plate	1.0	4.0 – 8.0 x 10 <sup>4</sup>
96-well Plate	0.3	1.2 – 2.4 x 10 <sup>4</sup>

## Procedure

### For suspend cells

The following protocol is given for transfection in 6-well plate. The protocol can be scaled up or down according to culture volume.

Cell Seeding: Suspension cells are typically seeded the day of the transfection at a density of 1.0~3.0 x 10<sup>5</sup> cells per ml of culture. For optimal transfection conditions with T-Pro Transfection Reagent II, seed the number of cells adapted to the culture vessel format according to Table 3.

- 1 For each well, dilute 5 µg of DNA into 50 µl of Serum-free Medium with High Glucose. Vortex gently and spin down briefly.
- 2 For each well, dilute 30 µl of T-Pro Transfection Reagent II reagent into 100 µl of Serum-free Medium with High Glucose. Vortex gently and spin down briefly.
- 3 Add the 130 µl T-Pro Transfection Reagent II solution immediately to the 55 µl DNA solution all at once.
- 4 Vortex- mix the solution immediately and spin down briefly to bring drops to the bottom of the tube.
- 5 Incubate for 15 minutes at room temperature.
- 6 Add the 200 µl T-Pro Transfection Reagent II / DNA mixture drop-wise onto the serum-containing medium in each well, homogenize the mixture by gently swirling the plate.
- 7 Incubate at 37 °C and 5% CO<sub>2</sub> in a humidified atmosphere.

Transfection experiments are usually stopped after 24 to 48 hours and gene activity assessed. Cells growing in suspension are collected by centrifugation at 800 x g and then resuspended in the desired medium or buffer.

**Table 2.** Recommended Amounts for Different Culture Vessel Formats.

<b>1</b>	<b>Subculture cells the day before transfection</b>				
		Adherent cells			Cells in suspension
	Culture dish diameter	35 mm	60 mm	100 mm	-
	Cell number (yielding 60-80% confluency)	1.0-3.0 × 10 <sup>5</sup> /dish	0.5-1.0 × 10 <sup>6</sup> /dish	1.5-2.5 × 10 <sup>6</sup> /dish	1.0 -3.0 × 10 <sup>5</sup> /ml
	Culture medium volume	1-2 ml	5-6 ml	10-14 ml	2-5 ml
<b>2</b>	<b>Prepare the T-Pro Transfection Reagent II/DNA mixture</b>				
2a	Dilute DNA to a final volume of	~2.5 µg 25 µl	~5 µg 50 µl	~75 µg 75 µl	~5 µg 50 µl
2b	Dilute T-Pro Transfection Reagent II to a final volume of	15 µl 50 µl	30 µl 100 µl	45 µl 140 µl	30 µl 100 µl
	Mix the DNA (2a) and T-Pro Transfection Reagent II (2b) solutions and incubate			<b>10-15 min</b>	<b>15-25°C</b>
<b>3</b>	<b>Transfect cells</b>				
	Replace culture medium by fresh medium Containing the T-Pro Transfection Reagent II /DNA mixture			<b>3~24 h</b>	<b>37°C</b>
<b>4</b>	<b>Further culture cells according to the experimental design</b>				<b>37°C</b>

**Table 3.** Recommended Number of Suspension Cells to Seed.

Culture Dish	Number of Cells
96-well plate	2 × 10 <sup>4</sup> - 5 × 10 <sup>4</sup>
48-well plate	5 × 10 <sup>4</sup> - 1 × 10 <sup>5</sup>
24-well plate	1 × 10 <sup>5</sup> - 2 × 10 <sup>5</sup>
6-well plate	2 × 10 <sup>5</sup> - 5 × 10 <sup>5</sup>
35 mm dish	5 × 10 <sup>5</sup> - 2 × 10 <sup>6</sup>
60 mm dish	2 × 10 <sup>6</sup> - 5 × 10 <sup>6</sup>
100 mm dish	5 × 10 <sup>6</sup> - 1 × 10 <sup>7</sup>